# **Sound Generation**

User Manual

Developed by: Camille Toarmino

For questions:

camilletoarmino@gmail.com


716-957-0746

# Contents:

- I. Creating a pure tone
- II. Creating broadband noise
- III. Creating a range of tones
- IV. Creating narrowband noise
- V. Loading .wav files
- VI. Saving data


### I. Creating a Pure Tone

- 1. Under CREATE, check the box that says 'Pure Tone'.
  - Make sure that there are no other boxes checked!
- 2. Enter a number in kHz for the frequency of the tone in the box that says 'Constant Freq' under PARAMETERS.
- 3. Enter the duration of the tone is milliseconds in the box titled 'Duration' under PARAMETERS.
- 4. Enter the number of iterations in the 'Repeats' box. If only presenting the stimulus one time, enter '1'.





#### Additions:

- To modulate the amplitude of tone, enter a number in Hz in the first box titled 'SAM'.
- To vary the amplitude modulation within a range, enter the low number in the first box and the high number in the second box.
  - Check the 'VARY' box next to SAM.
  - See more on amplitude modulation in section VI.
- If you have multiple iterations of the stimulus, you might want to add an ISI between each iteration.
  - o Enter a number in milliseconds in the ISI box.
- To randomize the presentation of multiple stimuli, check the 'Randomize' box next to the 'Generate' button.




## II. Creating Broadband Noise






### III. Creating a range of tones




#### \*When entering the range of tones:


- (1) The default is in octaves. If you enter a range that isn't naturally in octaves, it will only play to the last octave. If you choose force steps and enter a number, then the range will be split by the number of steps.
- (2) If you want to repeat an entire range of e.g. 4 tones 3 times, then multiply those numbers and enter 12 in the repeats box.

#### IV. Narrowband noise







## V. Loading .wav files



### VI. Saving Data



Every time you create a stimulus set and broadcast it, a file is automatically saved to MATLAB with the following cells: 1) start times, 2) stimulus ISI, 3) stimulus duration, and 4) stimulus type.



2.Option to save the entire log as a back up. Enter name of the .mat file. It will not be split up the way the automatically saved files are. It's sole purpose is to use as a reference and back up.

 Enter a filename for saving. Each file created will be noted with a number. To make it easier to find, I suggest making the filename something related to the stimulus set